

1-877-934-2336
www.ivoclarvivadent.com

Laboratory _____
 Contact _____
 Address _____
 Phone _____ email _____
 Case ID _____

Restoration Type:

- Crown Bridge
 Coping Framework

Required materials for submission:


- Stone or plastic models with opposing
- Pinned, cut and trimmed
- Articulated with plastic hinges or in plastic trays
- No die hardener or die spacer applied on the dies
- Do not send impressions

Material Type:

- IPS e.max® CAD
 Translucency (circle one): HT MT LT MO
 Zenostar® MT (crowns & bridges)
 Zenostar® T (crowns & bridges)
 Zenostar® MO (frameworks)
 Telio® CAD (temporary crowns & bridges)
 Shade Options: A1 A2 A3 A3.5 B1 BL3

Shade: _____

Notes:



Canada
 Ivoclar Vivadent, Inc.
 1 Riverside Drive West, Suite 402
 Windsor, Ontario N9A 5K3
 1-877-934-2336